





# MATTSON GENESIS FT-IR OPERATING INSTRUCTIONS

1. **Control Panel** should appear on the screen. If it doesn't, press the  icon. See diagram below.
2. Click **options** in the control panel.
  - a. click **LOAD METHOD**
  - b. click **tmp.ini** then **OK**
  - c. click **options**
3. With nothing in the sample compartment, click **Background** in the **Control Panel**.
  - a. press the green **scan** button in the **Control Panel**
  - b. click **OK** in the **Ready for Background Scan** box
  - c. click **OK** for title.



4. Place sample in the sample compartment.
  - a. click **scan** in the **Control Panel**
  - b. click **OK** in the **Ready for Sample Scan** window
  - c. type in what you want for a **title** when prompted and then **OK**
5. Data will appear, filling the entire window.
  - a. press **ctrl p** ---the peak picking window appears
  - b. click **OK**
  - c. a red box appears at the top of the computer screen instructing to **pick 2 points for threshold. DO NOT** press **OK**.
  - d. Put the pointer in the data window near the left Y axis at a height below which all the peaks will be picked. Pick a second point near the right side Y axis and click for second point.
  - e. The **Peak Report** window will appear. Click **X** in the upper right corner of the window to exit.
  - f. The peaks should all be labeled.
  - g. In the data area, right click the mouse twice. A new window, **DATA AREA** appears.
  - h. Change the value for **Bottom**, so all peak numbers show, i.e. -10. Click **OK**. You may have to repeat this step until all the peak numbers show.
  - i. Press **F7**, click **PLOT**, click **Done**.
  - j. Click the  icon in the tool bar to erase you spectrum.
  - k. Click the  icon so the instrument is ready for the next user.

10  
20  
40  
70  
100  
140  
180

6. **DO NOT** exit the **WinFIRST** window.

